

Using the DaRotary Program

The DaRotary data-collection program for the TI Graphing Calculators is used to collect, examine, analyze, and graph rotary motion sensor data. DaRotary is a single program, not a group.

**The DaRotary program for the TI-92 calculator, however is a group (DaRotary.92g). It contains the darotary.p and darotarz.p subprograms.*

Use the link cable to connect the LabPro or CBL 2 to the TI Graphing Calculator. Firmly press in the cable ends. Turn on the calculator. Follow these steps to start the DaRotary program on your calculator:

TI-73, TI-83 or TI-83 Plus Calculators

Press **PRGM** and then press the calculator key for the *number* that precedes DAROTARY (usually **1**). Press **ENTER** and wait for the main screen to load.

TI-86 Calculators

Press **PRGM** and then press **F1** to select <NAMES>, and press the menu key that represents DaRotary. (<DAROT> is usually **F1**). Press **ENTER**, and wait for the main screen to load.

TI-89, TI-92, or TI-92 Plus Calculators

Press **2nd** [VAR-LINK]. Use the cursor pad to scroll down to “DaRotary”, then press **ENTER**. Press **)** to complete the open parenthesis that follows “DaRotary” on the entry line and press **ENTER**. Wait for the main screen to load.

Main Screen

The main screen of DaRotary is shown at the right. The top half of the screen shows the current sensor setup and data collection mode. The portion below the double bar displays the menu options. The current status of the rotary motion sensor is shown. This status display is useful to verify that the sensor is functioning properly

```
DIG1:ROTARY(REV).00B
MODE:TIME GRAPH-30
-----
1:SETUP      4:ANALYZE
2:START     5:QUIT
3:GRAPH
```

SETUP

Choose SETUP from the main screen to select which rotary motion sensor you wish to use, add another sensor, and change the data-collection mode.

```
DIG1:ROTARY(REV)
CH 1:
MODE:TIME GRAPH-30
-----
1:OK
2:CALIBRATE
```

DIG 1

Here you can change various details about the rotary motion sensor. Using the up or down arrow keys on your calculator, position the arrow by the feature you want to change and press the **ENTER** key.

```
DIG 1:
TYPE:VERNIER/PASCO
RESOLUTION:LOW
UNITS:ROTARY(REV)
VEL/ACCEL:IGNORE
RADIUS(MM):27.5
-----
1:RETURN TO SETUP SCREEN
```

TYPE

Select which rotary motion sensor you are using: Vernier, PASCO, Team Labs, or CSMT

```
SELECT SENSOR TYPE
-----
1:VERNIER/PASCO
2:TEAM LABS (PSL)
3:CSMT
```

RESOLUTION

Press the **ENTER** key to toggle the resolution between low and high.

UNITS

Select the number corresponding to the unit of measurement you will be using.

```
SELECT SENSOR UNITS
-----
1:ROTARY(REV)
2:ROTARY(RAD)
3:ROTARY(DEG)
4:LINEAR(CMM)
5:LINEAR(CM)
```

VEL/ACCEL

Press the **ENTER** key to record the velocity and acceleration of your data.

RADIUS (MM)

Here you can change the radius of the pulley being used.

```
PULLEY RADIUS
IN MILLIMETERS:
?■
```

CH 1

Here you can select an additional sensor to be connected to Channel 1. After selecting the sensor, you will be returned to the Setup screen where you can calibrate it. This is especially helpful when selecting a custom sensor.

```
SELECT SENSOR
-----
1:TILIGHT
2:VERNIER LIGHT
3:DUAL RANGE FORCE
4:MAGNETIC FIELD
5:CUSTOM 0 TO 5V
6:CUSTOM -10 TO 10V
7:NONE
```

MODE

The default data collection mode for the rotary motion sensor is Time Graph (collecting 120 samples, one every .25 seconds). Position the cursor using the arrow keys opposite MODE and press **ENTER** to see the Select Mode screen. Press the number corresponding to the mode you want.

```
SELECT MODE
-----
1:TIME GRAPH
2:EVENTS WITH ENTRY
3:SELECTED EVENTS
4:RETURN TO SETUP SCREEN
```

TIME GRAPH

Choose CHANGE TIME SETTINGS to enter a new time interval or number of samples. Choose OK to return to the Main Menu.

TIME GRAPH SETTINGS	
TIME INTERVAL:	.25
NUMBER OF SAMPLES:	120
EXPERIMENT LENGTH:	30
1:OK	
2:CHANGE TIME SETTINGS	

EVENTS WITH ENTRY

Data collection occurs one point at a time and only when the button on the calculator is pressed. You are then prompted to enter the corresponding value.

SELECTED EVENTS

This mode collects data from each channel in response to pressing the key, but the user does not enter any values. Point numbers are plotted vs. data.

START

Select Start to begin collecting data using the existing sensor setup and data collection mode.

GRAPH

The GRAPH screen automatically appears when you are finished collecting data or when GRAPH is chosen from the main screen. To view a graph, select the appropriate channel from the list by using the and keys on the calculator and pressing .


DIG-ROTARY(REV)	
1:RETURN TO MAIN SCREEN	
2:SELECT REGION	

If you want to view another graph, press . The graph screen appears again so you can choose another channel.


SELECT REGION

This option is used to remove unwanted data from lists, *i.e.* data that are outside the region of interest. This option may be necessary for data analysis. Data outside the selected region are permanently deleted.


- Select the left bound, as viewed on the bottom, using the and keys. When the values you want for the left bound is displayed, press .



- Repeat for the right bound and press .



- The new graph will automatically be displayed.



ANALYZE

The ANALYZE option on the main screen produces this list of options. You can perform a curve fit, determine statistics over a region of data, and determine the integral over a region of data.

ANALYZE MENU	
1:	RETURN TO MAIN SCREEN
2:	CURVE FIT
3:	STATISTICS
4:	INTEGRAL

CURVE FIT

Choosing CURVE FIT from the ANALYZE OPTIONS menu produces a list of curve fit options.

SELECT CURVE FIT	
1:	LINEAR(ROTARY VS TIME)
2:	LINEAR(VEL VS TIME)
3:	QUAD(ROTARY VS TIME)
4:	POWER(CH1 VS ROTARY)
5:	RETURN TO ANALYZE OPT.

After a curve fit is performed, the fitted curve is drawn as a continuous line, and the data points are shown with box point protectors.

STATISTICS

This option finds the mean, min, max, standard deviation, and number of data points of a selected region. After selecting this option, choose a graph from the displayed list. Use the arrow keys to move the cursor and select the left bound of the region and press **ENTER**. A vertical line will be drawn on the graph. Now move the cursor to select the right bound and press **ENTER**. After the calculator has finished determining the statistics, a message will be displayed instructing you to press **ENTER** to continue. The numerical results will be displayed. Calculating statistics does not delete data.

SELECT GRAPH	
1:	ROTARY
2:	VELOCITY
3:	ACCELERATION
4:	CH1
5:	RETURN TO ANALYZE OPT.

INTEGRAL

The integral function is used to integrate a section of the graph. After selecting this option, choose a graph from the list. Use the arrow keys to move the cursor and select the left bound of the region and press **ENTER**. A vertical line will be drawn on the graph. Use the cursor to select the right bound and press **ENTER**. The numerical results will be presented.

SELECT GRAPH	
1:	ROTARY
2:	VELOCITY
3:	ACCELERATION
4:	CH1
5:	RETURN TO ANALYZE OPT.

QUIT

Select Quit to leave DataMate and return to the home screen of the calculator. On the TI-89, TI-92 and TI-92 Plus calculators you must press F5 or HOME to return to the home screen. To see the last graph select zoomstat as appropriate for the calculator model.