

Logger *Pro* 3 Release Notes

Version 3.8.4

February 1, 2011

Contents

System Requirements

Install Notes

What's Changed in Logger *Pro* 3.8.4

Logger *Pro* 3.8.4 Known Issues

File Compatibility

What Changed in Earlier Logger *Pro* 3 Releases

What Changed from Logger *Pro* 2 to Logger *Pro* 3

LabPro Operating System Changes

LabQuest Operating System Changes

General System Requirements

Basic data collection or graphing, with no video capture

- Windows XP SP2/SP3 (32 bit) with Pentium processor, 500 MHz or better, 256 MB RAM, available USB or serial port. Windows Vista or Windows 7 (32 or 64 bit) with OS compatible hardware. Apple's QuickTime version 7.1.6 or newer required for video playback.
- Macintosh (PPC & Intel) OS X 10.4.1, 10.5.3, 10.6, or newer with available USB port. We recommend that Mac users use 10.5.8 or better and upgrade their systems through Apple's software update on a regular basis for the optimal experience in using Logger *Pro*.
- Logger *Pro* requires approximately 380 MB. About 150 MB additional temporary disk space is needed during the installation process; you should have at least 830 MB free space prior to installation to allow the installation to complete successfully and retain some space for user files. Video playback requires QuickTime – another 410 MB should be reserved prior to a Windows installation.
- Computer display resolution 1024x768 or higher is strongly recommended. Lower resolutions may result in important controls being hidden
- For video playback or extended data collection a faster processor and more memory are suggested, such as a 2 GHz (867Mhz G4 for Mac) processor with 2 GB RAM. QuickTime 7 is required for video playback.
- Data collection requires one or more of the following interfaces: LabQuest Mini, LabQuest, LabPro, Go! Temp, Go! Link, Go! Motion, WDSS, Vernier Spectrometers, select Ocean Optics Spectrometers, select Garmin GPS receivers, Vernier GPS, Mini GC, Watts Up? Pro, Ohaus Scout Pro or Adventurer Pro balance.
- Previous releases of Logger *Pro* are available for special applications. Contact Vernier for more information.
 - Video resolutions less than 1024x768 are fully supported by Logger *Pro* 3.4.6 for non-video applications.
 - Window 2000 is supported by Logger *Pro* 3.6.0.

- Windows 98 SE & Me are supported by Logger *Pro* 3.4.6.
- Mac OS 9 is supported by Logger *Pro* 3.4.2; OS X 10.2 is supported by Logger *Pro* 3.4.5; OS X 10.3 is supported by Logger *Pro* 3.6.1
- Windows NT is supported by earlier versions of Logger *Pro* 3. Windows NT does not support any USB devices.
- Windows 95 and Mac OS 8 are supported by Logger *Pro* 2.2.1.

Video Capture Functions

- Windows XP, Vista, or Windows 7 with DirectX 9 and OS compatible drivers installed.
- Macintosh OS X 10.4.1 or better. OS X 10.4.8 or newer is required for Proscope HR use on Intel Macs.
- For Windows: 2 GHz single-core processor or any multi-core processor and 1GB RAM, and at least 500 MB available hard disk space for temporary storage.
- For Mac: 867 MHz G4 processor and 1 GB RAM, and at least 500 MB available hard disk space for temporary storage.
- Routine video capture and storage will require additional disk space over time.
- Screen resolution of at least 1024x768.
- Compatible camera.
- A minimal system may sometimes exhibit skipped or duplicated frames in video capture. Faster systems will exhibit fewer such problems. For most purposes these duplicated or skipped frames will not be visible.

Wireless / Bluetooth with Vernier Wireless Dynamics Sensor System (WDSS)

- For Windows, WDSS requires Windows XP SP 2 or newer, Vista, or Windows 7 with the Microsoft Bluetooth Driver and a Bluetooth radio compatible with the installed driver (built-in or add-on radio). Certain “Widcomm” or “Broadcom” Bluetooth radios are supported. Vernier recommends using the Microsoft Bluetooth driver whenever possible for such radios; however, latest Broadcom drivers and radios have been shown to work with Logger *Pro*.
- For Macintosh, WDSS requires OS X 10.4.1 or newer (10.4.11 or newer recommended); a Mac-compatible Bluetooth radio (built-in or add-on).
- For recommended Bluetooth radios or other Bluetooth/WDSS related information, visit www.vernier.com/bthelp

Install Notes

There is an interaction between the Logger *Pro* installer and Texas Instruments TI Connect uninstaller. If TI Connect has been installed and subsequently removed, Logger *Pro* installer will refuse to run. To install Logger *Pro* 3 in this situation, first reinstall TI Connect 1.6 or newer. You can download an installer from <http://education.ti.com>.

Logger *Pro* 2, Logger Lite, Graphical Analysis, LabQuest Emulator, and Logger *Pro* 3.8.4 will coexist on a single computer.

Logger *Pro* is supported in the following languages: Arabic, Czech, English, German, Spanish, French, Italian, Dutch, Polish, Portuguese, Russian, Simplified and Traditional Chinese, and Korean.

Video capture and picture analysis functions in Logger *Pro* require QuickTime 7 installation. Normally a Logger *Pro* installation includes QuickTime 7 on Windows computers (Mac computers' QuickTime installations are left alone by Logger *Pro* installers). However if the Silent Installer method is used or if the person installing Logger *Pro*

opts to omit QuickTime installation, you will need to install QuickTime 7 to use the picture and video functions in *Logger Pro*.

Vernier's *Logger Pro* installers for Windows provides an Apple installer for QuickTime as a courtesy to the users of *Logger Pro* in support of *Logger Pro*'s video functionality. The specific QuickTime version provided by Vernier's Windows installation process was current as of the date that version of *Logger Pro* went to Beta release status. After the start of a *Logger Pro* Beta, Apple may independently elect to release a new Windows version of QuickTime. Thus it is possible that when *Logger Pro* is installed one of five conditions may arise regarding QuickTime.

- No version of QuickTime is installed on the target Windows machine. In this case, the user installing *Logger Pro* is encouraged to allow the installation of QuickTime by the *Logger Pro* installation sequence.
- An older version of QuickTime already exists on the target Windows machine. In this case, please allow the QuickTime installer included with *Logger Pro* to proceed with the installation of QuickTime.
- The same version of QuickTime is already on the target Windows machine. The *Logger Pro* installer detects this and skips the QuickTime installation. The *Logger Pro* installer will complete normally.
- A new version of QuickTime is already on the target Windows machine. The *Logger Pro* installer detects this situation and skips the QuickTime installation. The *Logger Pro* installer will complete normally.
- After *Logger Pro* installs a version of QuickTime, Apple may detect later that the version is older than an available new version. This is not a problem with Vernier's installation sequence but is normal and expected if Apple releases a new version of QuickTime after *Logger Pro*'s installer was created. Users are encouraged to update as Apple indicates. This is increasingly likely if there has been a delay in installing a new version of *Logger Pro* after it is first made available; but, it may occur even immediately upon the official release of *Logger Pro* since *Logger Pro* installers are not updated for a few weeks prior to release to ensure testing stability for *Logger Pro*'s video functionality.

In the rare event that an Apple update to QuickTime adversely affects *Logger Pro* functionality, please contact Vernier Technical Support.

Windows Network Installation/Remote Deployment

- An installer for Windows that does not prompt for user input is available by request. This is useful for remote or network managed installations of *Logger Pro*, as this silent installer can be caused to run on a remote computer, and the installation will proceed. Installation on Vista and Windows 7 systems may require User Access Control (UAC) changes (temporarily) in order to complete the installation via the silent installer. Such UAC changes are ultimately the responsibility of the computer's administrator and/or network management tools. The silent installer does not include QuickTime. If you require picture and video playback on your remote or network managed Windows installation, you must install QuickTime as a separate operation. Write to info@vernier.com for more information.
- We do not recommend installing *Logger Pro* to a network location for shared access, because the necessary drivers will not be present on client machines.
- We do not recommend re-packaging a *Logger Pro* installation for remote deployment, because driver installation requires that a secondary installer execute at install time. Please use a Vernier supplied silent installer with your remote deployment tools to accomplish needed remote installation scenarios.

What's Changed in Logger *Pro* 3.8.4

New Device/Interface/Sensor Support

- Improved WDSS Bluetooth Support – newer Broadcom (Widcomm) Bluetooth radios support.
- Supports new anemometer sensor
- Supports new Melt Station
- Supports new 30V sensor

New Features

Data Markers enable annotations of data events during a time-based data collection. The Data Markers are associated with the actual data cell(s) taken at the user selected time during the collection. Such markers might for instance be associated with a change of phase during a melting session with the new Melt Station or any other sensor where you may wish to note an event during the collection. Hit the <TAB> key during a data collection places a data marker and displays an associated helper box for further use. Double-clicking the associated helper box allows you to enter text associated with that mark.

Data Tagging enables annotations of specific cells or ranges of cells after a collection or for collections utilizing non-time-based modes. Similar to Data Markers, Data Tagging can associate additional information with specific data elements.

Note that Data Markers and Tags are distinct from text annotations in that the marker or tag is associated with a distinct cell in the data table and not just a location on a graph as is the case for text annotations. Use Data Markers as part of a time-based collection during the collection. Use Data Tagging after a collection. Data Tagging can be modified repeatedly as needed and re-associated with different data cells. Data Markers are restricted to the original data cell or cells with which they are first associated during a collection and may not be re-associated with other cells.

More video analysis options are available. Video frames may be skipped while adding points. The first point marked may set time zero. Center of mass studies can be done with video analysis.

Other Changes

The Logger *Pro* Windows installer has been rewritten to better support international or second language installations, enable better Windows support of Broadcom based Bluetooth radios, enable better installer behavior when the same or newer version of QuickTime is already installed, rectify some driver issues that plague some XP installations, and provide MSI based “silent installation” for improved remote/network managed deployment usage.

Czech, Dutch, and Polish were added as supported languages.

Logger Pro 3.8.4 Known Issues

Disabling Auto Behavior

Logger Pro 3.8.4 can be configured to not automatically setup for collection based on the sensor connected in support of certain curricula; however, many devices do not yet allow “manual” setups. Examples include Spectrometers, Ohaus balances, Go!Temp, and Go!Motion. Only sensors which use the LabPro, LabQuest, LabQuest Mini, or Go!Link may be used in the manual setups.

LabQuest file Compatibility

Logger Pro 3.8.4 can open most LabQuest files correctly. Data columns will be translated correctly. However, in some cases there will be errors in calculated column translation, breaking the data import. In most cases you can manually recreate the calculated columns in Logger Pro to recover the information or use fully supported equation forms in files meant for interchange between Logger Pro and LabQuest.

Similarly, LabQuest App can open many Logger Pro files when transferred to the LabQuest using the LabQuest Browser in Logger Pro. Elements not supported on LabQuest will be ignored, such as text annotations or movies. While LabQuest App does support calculated columns, the possible calculations are much more limited than in Logger Pro. As a result, some calculated columns will not be translated for use in LabQuest App.

LabQuest Device Compatibility

Logger Pro 3.8.4 can be used with LabQuest 1.0, 1.1, 1.1.1, 1.2, 1.3, 1.4, and 1.5 software. If a LabQuest with 1.0 software is detected, a warning will be displayed advising that the LabQuest unit be upgraded. This message will be displayed every time the LabQuest is detected, but the LabQuest can still be used for most applications. If a Logger Pro feature requiring newer LabQuest software is used, unexpected results such as freezing, crashing, or data loss may occur. Power Amp and Function Generator functionality is only available on version 1.3 of LabQuest or newer. We strongly recommend updating LabQuest to 1.5 or newer. See www.vernier.com/labquest/updates/ for information on updates.

Windows 7

The behavior of USB devices after a computer sleep/wake cycle on Windows 7 may differ from prior Windows operating systems. For example an interface connected to Logger Pro (not actively collecting data) may lose the connection to the interface when the computer wakes from a sleep state.

Typically the computer enters a sleep state as the result of the computer's energy settings. When Logger Pro is actively collecting data the program prevents the computer from sleeping; however, if the interface is merely connected but not collecting Logger Pro will allow the sleep cycle to start.

If your Windows 7 installation was installed fresh (not an upgrade from Vista) then you will likely see the interface fail to reconnect upon computer wake. Upgrades from Vista retain the Vista behavior and the interface will reconnect automatically upon computer wake. If you encounter problems with USB connected interfaces in Windows 7 upon a computer waking from sleep, it is recommended that you disallow sleep, disconnect interfaces

when done with active collection, or set the time to commence sleep well beyond the anticipated use of Logger Pro on that computer.

In the event an interface fails to reconnect, shut down Logger Pro, disconnect the interface, reconnect, and invoke Logger Pro again. The interface should then connect without a problem.

Vista & Windows 7

Logitech 4000 cameras are not supported; use Logitech 5000 model instead.

Uninstalling 3.8.4 from the control panel will complain about an unknown publisher. This does not affect upgrades from prior versions. To work around the issue, simply ignore the unknown publisher warning and allow the uninstaller to remove Logger Pro 3.8.4 as desired.

Windows XP

Uninstalling 3.8.3 will display a dialog during the uninstalling of experiments that will be unreadable due to missing fonts unless extended language support has been installed on the machine. This is a InstallShield and Windows XP interaction; the uninstall will complete normally and no user action is required.

Windows 2000

Windows 2000 is supported by Logger *Pro* 3.6.0; please use 3.6.0 for use on Windows 2000. The LabQuest interface is not supported on Windows 2000 in any version of Logger *Pro*.

Intel Macs

Original Proscopes on Intel Macs requires use of Rosetta emulation mode and can be problematic; Proscope HR is recommended for Intel Macs.

Proscope usage on Intel Macs requires using Mac OS X 10.4.8 or newer. If the Apple built-in iSight camera is present it may initially appear to preclude using the Proscope models. To use an external camera on Intel Macs running 10.4.9, launch Logger *Pro* and insert a video capture object. Logger *Pro* will use the iSight. Click the Options button, and then the Camera Settings button. Click the Source tab. There you can choose to use the Proscope or other external camera as long as the drivers for that camera are present. Details may be different for other versions of the OS. Please contact Vernier if you encounter difficulties and reference TIL #1643.

Spectrometers

Although multiple spectrometers may be connected at the same time and used with Logger *Pro*, such use is not supported by Vernier. We recommend that only one spectrometer be connected at a time.

You may need to reboot on Vista after connecting your Ocean Optics spectrometer for the first time. This is normal.

If you experience connection issues with the USB2000 or other Ocean Optics spectrometers after installing the latest *Logger Pro* and rebooting the machine (particularly on Vista), please run Windows Update to ensure that Vista is fully updated with all Microsoft's updates for the operating system and try connecting after updates are complete. Vista 32 bit systems that are not up to date with Microsoft's updates may exhibit connection problems particularly with the USB2000 unit.

GPS

At this time only Garmin devices are supported with the Garmin data transfer protocol. NMEA data protocols are not supported in 3.8.4. Only select Garmin units and the Vernier GPS devices are supported at this time; please see www.vernier.com/gps for supported models.

COM interfaces are shown for connected VGPS devices. Attempting to re-connect to a Vernier GPS (VGPS) manually via a COM port (after disconnection for whatever reason) will not yield a usable connection and may lock up or crash *Logger Pro*. If there is a need to reconnect to the VGPS it is currently required for you to use the Vernier GPS port designation in the Experiment>Connect Interface>GPS> menu.

VGPS connection to Mac 10.6 (Snow Leopard) mistakenly will display a networking dialog. Users of VGPS on Snow Leopard may simply cancel from that dialog and use the device normally in *Logger Pro*.

Adventurer Pro

Pressing the "Unit" button on the Ohaus scale may lead to incorrect values being reported by *Logger Pro*. Use the Calibrate dialog in *Logger Pro* to set the units.

Watts Up? Pro

Although the device contains on-board memory and internally logs data values, *Logger Pro* does not support retrieving that stored data (remote data retrieval). To use the "Watts Up? Pro" device with *Logger Pro* you must collect data with the device connected.

Mac OS 10.3

Mac users of OS X 10.3 should use *Logger Pro* 3.6.1. You must use QuickTime 7.1.6 or newer. Specifically, QuickTime 7.1.5 and *Logger Pro* 3.6.1 are not compatible due to Apple specific issues on OS X 10.3.

General Issues and Suggestions for Logger Pro 3.8.4

- Currently you cannot insert the same movie in an experiment file (cmbl) more than once and have them share the same video analysis data set. Duplicate movies will have independent data sets.
- A known Apple QuickTime issue (versions prior to 7.5) causes cropping of videos in *Logger Pro* if the *Logger Pro* window is sized less than full screen. Video analysis will not be possible if the *Logger Pro* window is sized to less than full screen using obsolete QuickTime versions. Use full screen *Logger Pro* with

using videos unless you have upgraded your QuickTime installation to version 7.5 or greater on Mac OS X 10.4 or greater (10.3 with QuickTime 7.5 continues to exhibit the cropping issue).

- Calculated columns' equations based on two or more data columns measuring the same physical quantity (such as temperature) from sensor data collected using the LabQuest interface will not be preserved if saved as QMBL format and read by *Logger Pro*. For example, if you measure two temperatures with the LabQuest (standalone) and save a QMBL formatted file set with a calculated column which depends on those two temperatures, save, and reopen that QMBL formatted file in *Logger Pro*, the calculated column will not work as expected. The workaround is to re-create the calculated column in *Logger Pro*.
- Do not connect the USB cable to LabPro until it is powered and has given its start-up beeps.
- Help doesn't display correctly in Mac OS 10.4. This is a Mac OS defect; to display help properly update to 10.4.1 or newer.
- Meters monitoring Drop Counter sensors will display the last point taken.
- System response time may increase when collecting a large number of data points.
- When using a Stainless Steel Temperature Probe remotely, the default calibration will be used. A custom calibration can be performed, but will be used only when the collection is done while attached to a computer.
- When updating the LabPro OS on multiple LabPros, do not disconnect a LabPro until after you click OK on the dialog box announcing a successful update.
- User parameters cannot start with z.
- Photogate data collection can fail in digital events mode with a bad time value. This issue also affects earlier releases, and will be corrected in a LabPro firmware update. Contact Vernier Software & Technology at info@vernier.com for more information.
- On Macintosh only, picture analysis may not work if the inserted image is a .pct file. To avoid the problem, convert to .jpg in an image editing program.

Logger Pro File Compatibility

Logger Pro 3.8.4 can open files created in *Logger Pro* 3.0, 3.1, 3.2, 3.2.1, 3.3, 3.4, 3.4.1, 3.4.2, 3.4.5, 3.4.6, 3.5.0, 3.6.0, 3.6.1, 3.7, 3.8, 3.8.2, and 3.8.3. *Logger Pro* 3.8.4 can also open files created in *Logger Lite* (any version) and *Graphical Analysis* 3.0, 3.1, 3.2, and 3.4.

Files created by *Logger Pro* 3.8.4 can be opened in version 3.3 through 3.8.3, but if 3.8.4 features are used in the file, it will not be fully functional.

Files created by *Logger Pro* 3.8.4 cannot be opened in *Logger Pro* 3.0, 3.1 or 3.2.1. *Graphical Analysis* cannot open any *Logger Pro* files.

Specific File Issues

- Files created in *Logger Pro* 3.8.4 containing calculated columns with the degree symbol or Greek micro (μ) in the column name will not open correctly in *Logger Pro* 3.3. To avoid the problem, remove the degree symbol or Greek micro character while in version 3.6.2. Or, in version 3.3, find the $\hat{\text{A}}^\circ$ or $\hat{\text{A}}\mu$ you see in the column name and remove the $\hat{\text{A}}$.
- Files created in *Logger Pro* 3.3 or earlier containing column names that use the pipe character (|) will not open correctly in *Logger Pro* 3.8.4. To resolve this problem, remove the pipe character from the column names while using *Logger Pro* 3.3 or earlier.
- Files created in 3.8.4 that include calculated columns referring to a specific data set will not open properly in *Logger Pro* 3.3.
- *Logger Pro* 3.8.4 continues to support auto-save files as a safeguard against unexpected interruption of data collection. Increasingly, experiment files can contain large numbers of samples that take a long time to save

to disk. To remedy this in 3.6.0 or later, auto-save files are now saved in a new internal format allowing, for example, files that formerly took 40 seconds to save, take only one second to auto-save.

Changes in Earlier Logger *Pro* Releases

Logger *Pro* 3.8.3

New Device/Interface/Sensor Support

- Latest VGPS devices (MTK) Support
- Ocean Optics Spectrometers supported on Windows 7 and Vista 64-bit operating systems.
- Drop Counter Automatic Configuration Feature Added

Other Changes

For the first time both English and international versions of Logger Pro 3 are concurrently released. English users in most cases will not notice any difference. However, if an English user wishes to enable one or more available alternate languages such as Spanish or French along with English it is relatively simple to do so by creating additional desktop shortcuts that specify a “/Lang” parameter to the Logger Pro 3 executable. For example, Spanish Logger Pro 3 on an English system may be enabled via a copy of the regular desktop shortcut and adding “/Lang=es” to the Target line in the shortcut. – contact Vernier for more information.

Current international users will enjoy more timely and seamless support of their native language on their choice of Windows operating systems (International Logger *Pro* 3 is not supported on Mac OS X) International users wishing to also access the English version (or other installed languages) of Logger *Pro* 3 on their native systems may do so by creating an additional shortcut(s) as above (e.g. “/Lang=en” for English). Logger *Pro* 3.8.3 replaces all prior versions of Logger *Pro* International.

Graphing optimized for user entered data when no interface is connected; no connecting lines are drawn and the scale is set to 0 to 10.

Logger *Pro* 3.8.2

New Device/Interface/Sensor Support

- Vernier LabQuest Mini (see the Vernier website for more information: <http://www.vernier.com/lqmini>)
- SpectroVis Plus
- Windows 7 64-bit operating systems support.
- Mac OS X 10.6 (Snow Leopard) support.

Other Changes

The location for Experiment templates and examples is changed for 3.8.2 to mesh with operating systems providing a “Common Application Data” location and avoid problematic automatic Windows behavior like “Shadow Copies” resulting from Logger Pro writing to “Program Files” routinely outside of the installation process. As a result, first time opens will default to the Common Application Data location (which varies by Operating System) for access to the provided Experiments. New users of Logger Pro will not notice this change; however, long time users of Logger Pro will no longer find the “Experiments” folder in the “Program Files” or “Applications” folders. The new locations are:

- /Library/Application Support/Logger Pro for Mac OS X
- [Startup Drive Letter]:\Documents and Settings\All Users\Application Data\Vernier\Logger Pro 3\{language code} for Windows XP
- [Startup Drive Letter]:\ProgramData\Vernier\Logger Pro 3\{language code} for Vista & Windows 7 systems

For the convenience of long time users of Logger Pro 3, an alias or shortcut that points to the “Experiments” folder may still be found in the “Program Files” or “Applications” folder.

Logger Pro 3.8

New Device/Interface/Sensor Support

- Vernier Mini GC (see the Vernier website <http://www.vernier.com/probes/gc-mini.html> for more information)
- LabQuest Function Generator (with LabQuest 1.3)
- LabQuest Power Amplifier (with LabQuest 1.3)
- Vernier VGPS-PL
- Watts Up? Pro

Windows 7 & Mac OS X 10.6 OS Provisional Support

Logger Pro is now provisionally supported on Windows 7 32-bit and 64-bit OS and Mac OS X 10.6 (Snow Leopard) Note that 10.6 is only available on Intel architecture Macs. Logger Pro 3.8 was tested on pre-release versions of both Windows 7 and Mac OS X 10.6 and found to be compatible. It is possible however that final released versions of these operating systems from Microsoft and Apple could introduce changes from the pre-released versions tested during Logger Pro 3.8 development and evaluation that could be detrimental to Logger Pro 3.8's operation. Windows 7 and Mac OS X 10.6 users are cautioned to perform their own evaluation of Logger Pro 3.8 on the released versions of these operating systems prior to committing to any large scale deployments of Logger Pro 3.8 on these new operating systems.

New Features

- Peak Area Measurement (Peak Integration) When you choose Peak Integration from the Analyze menu, a dialog will be displayed. Use the graph in the Peak Integration dialog to select the left and right boundaries of the peak you want to integrate. You can auto-scale, zoom in, and zoom out using the buttons. This feature primarily supports the analysis of chromatographs from the Vernier Mini Gas Chromatograph (GC)
- Function Generation with connected 1.3 LabQuest. To set up the Function Generator, open the LabQuest sensor dialog by clicking the LabQuest icon in the toolbar or choosing Set Up Sensors > LabQuest from the Experiment menu.

- Power Amplification with connected 1.3 LabQuest and Vernier's PAMP. To set up the Power Amplifier, connect PAMP to LabQuest and then open the LabQuest sensor dialog by clicking the LabQuest icon in the toolbar or choosing Set Up Sensors > LabQuest from the Experiment menu.

Significant Defects & Issues Resolved

In prior releases, unplugging devices would not remove them from the Live Readout bar under some conditions. This could confuse some users that did not always carefully disconnect and connect devices at the beginning and end of Logger Pro sessions. Device signatures from the low level API were corrected to match what Logger Pro expected in these cases so that the confusing behavior no longer occurs.

In 3.7, VGPS users of Logger Pro in the eastern hemisphere could encounter longitudes that were incorrectly signed. Such users would have to introduce calculated columns to correct the sign prior to using such data. Logger Pro 3.8 corrects this error.

Using a heat pulser (HP-BTD) with the 1.2 LabQuest interface attached to Logger Pro 3.7 did not work. Using a 1.3 LabQuest interface with Logger Pro 3.8 corrects issues with the heat pulser via the LabQuest interface. LabPro users are unaffected.

LabPro users attempting to upgrade LabPro firmware with Logger *Pro* 3.7 “over the wire” firmware updates would fail with an assert. The issue is corrected in Logger *Pro* 3.8. Note that upgrades are enforced for LabPro firmware older than 6.26.

Logger *Pro* 3.7

New Device/Interface/Sensor Support

- LabQuest 1.2 devices are fully supported including sensor support for wide range temperature, heat pulser, and DCU
- Vernier GPS Device

64 bit Vista OS Support

Qualified and signed 64 bit Device drivers for LabQuest and LabPro interfaces are now available. The Logger *Pro* installer automatically installs these new drivers enabling full support of Logger *Pro* functionality on this 64 bit operating system.

Better Date & Time Display

Date and Time display on graphs is improved, providing better user feedback and the ability to plot time-based data in much more meaningful and useful ways.

Improved GPS Functionality

- Signal strength indication for VGPS devices in Live Readouts.
- Ability to select coordinate units (decimal degrees, degrees + minutes, UTM)

- Better automatic data columns based on collection mode / sensor configuration.

Changes in Auto Behavior

- An entry in the usersettings.txt file, “IDSX_MANUAL_SENSOR_CONFIG” may be used to specify if automatic setup of sensors should be turned off to support curricula that have students perform setup for data collection, etc.
- The usersettings.txt is found in the Vernier installation folder; you may find this in typical installations under your Local Disk C: > Program Files > Vernier Software > Logger Pro 3 > en-US > Support > Usersettings.txt.
- In the Usersettings.txt file find the IDSX_MANUAL_SENSOR_CONFIG entry and change the value to “1” to enable the manual mode.

Installer Packaging Improvements

The Windows installer has been migrated to use a MSI package which supports easier deployment by IT personnel to multiple computers by scripting. Contact Vernier for more information.

Significant Defects & Issues Resolved

- Fixed: Multiplication in LabQuest calculated column fails to transfer to LP
- Fixed: Leaving LabQuest Interface window displayed in Logger *Pro* will cause interface disconnects.
- Fixed: cannot correctly zero a motion detector when using LabPro

Other Issues Resolved in 3.7

- Fuller set of calculated column equations supported in Logger *Pro* / LabQuest data exchanges.
- Improved LabQuest remote setup.

Logger *Pro* 3.6.1

New Device/Interface/Sensor Support

- LabQuest 1.1 devices fully supported including sensor support for radiation and rotary motion
- Support for LabQuest 1.1 internal microphone

New Icon

Logger *Pro* 3.6.1 features new icons for the application and documents. The new icons are related to new icons used in Logger Lite 1.4.

Significant Defects & Issues Resolved

- Crashes with WDSS when used with WIDCOMM Bluetooth driver

Other Issues Resolved in 3.6.1

- Clarified wording in Remote Data Retrieval dialogs

Logger Pro 3.6.0

New Device/Interface Support

- LabQuest devices supported (Note: International users (3.6.0i) may encounter specific file interchange issues, see LabQuest 1.2 or 1.3 notes at the end of this document.)
- New SpectroVis Spectrometer supported
- Global Positioning System (GPS) devices supported on Mac.

New Program Features

- Auto-saving during collecting is significantly faster.

Significant Defects & Issues Resolved

- Improved reliability during sensor confirmation and LabPro firmware updates.
- “Rough and slow graphing” on Macs displaying pixel ratio resolutions other than 1:1.

Other Issues Resolved in 3.6.0

- On Windows, the installers now install a version (7.2) of QuickTime supporting Vista computers.

Logger Pro 3.5.0

New Computer Operating System Support

- Windows Vista is supported on Vista compatible hardware. LabPro drivers are signed for Vista.

New Device/Interface Support

- Logitech 5000 camera for video capture.
- Garmin GPS receiver models GPSmap 60CS and GPSmap 76S.
- New Vernier CO₂ sensors with selectable range.
- Vernier and Ocean Optics Spectrometers are now supported on Intel Macs.

New Program Features

- Improved axis based zoom in/out on graphs; uses and extends the “stretchy axis” concept.
- Mac allows customization of the toolbar to suit user preferences.
- Mac User Interface conforms much more closely to Apple Human Interface Guidelines and thus should be more similar to other Mac applications. Mac Logger *Pro* text handling now has access to more of Mac native text handling and features.
- Numerous spectrometer connection and functional improvements. Connection/disconnection is faster and connection does not require a separate scan operation. Peak resolutions are generally cleaner and noise is diminished. Calibrations may reliably take advantage of averaging.
- Automatic graph match exercise creation feature. Access in Settings for (file-name) in the File menu.
- Empty Latest data sets are now hidden on remote data retrieval.
- Movies can now be swapped without altering video analysis environment.
- Global, as opposed to per-user, preferences.
- Sensor confirmation dialog now allows connection of WDSS.
- Additional picture and movie file types available in Mac insert picture/movie object open dialogs.

Significant Defects & Issues Resolved

- Macintosh long-term data collection spanning more than 3 days is more reliable.
- Curve fits no longer require sorting data.
- Numerous graphics and print issues.
- Spectrometer calibration more robust.
- Improved stability during drop counter calibration.
- Improved reliability during sensor confirmation and LabPro firmware updates.

Other Issues Resolved in 3.5.0

- Radiation monitor and events with entry or selected events again works.
- Removed some limitations on columns choice for horizontal axis.
- Allow simultaneous heart rate monitor and blood pressure measurement.
- Restored robust movie file location search feature on Macs.
- User parameter units no longer lost.
- Data Pro Hot-sync from Palm on Mac again automatically retrieves data into Logger *Pro*.

Logger *Pro* 3.4.6

New Computer Hardware Support

- Intel Architecture Macintosh computers are supported natively for the first time in Logger *Pro* 3.4.6 through the new Apple Universal Binary.

New Device/Interface Support

- IR Temperature Probe
- Garmin GPS USB interface initially supporting GPS Garmin Models Vista Cx and GPSmap 76S

- WDSS WIDCOMM interface (alternate Bluetooth driver) for increased Bluetooth configuration options on Windows. Note: Vernier strongly recommends using Microsoft's Bluetooth drivers on Windows XP SP2 whenever possible.

New Program Features

- Time Offset option for curve fit
- Import Logger *Pro* experiment files; combine several students' results in one Logger *Pro* file.

Significant Defects & Issues Resolved

- Keyboard Triggers now work properly
- WDSS repeat remote setup no longer loses trigger setting
- GO!Motion dialog for keeping remote data on device is no longer missing
- Eliminated assert on video analysis WUtils.cpp when user deletes the time column
- Improved tracking of multiple objects in video analysis

Logger *Pro* 3.4.5

New Device Support

- Spectrometers powered by Ocean Optics (XP, Windows 2000, Mac OS X 10.3.9 or newer)
- Wireless Dynamics Sensor System (WDSS) (XP SP 2, Mac OS X 10.3.9 or newer)
- Global Positioning Satellite (GPS) Receivers by Garmin (Windows only)
- Adventurer *Pro*
- Soil Moisture Sensor
- Digital Video compatible cameras and devices.

New Program Features

- Log Graphs (Single and Double Y Axis)
- Point and Click video / data synchronization
- Full spectrum mode, and visible spectrum graphing with Spectrometers
- Data synchronization with multiple WDSS devices (facilitates dynamics experiments that were problematic in the past)
- Export Latitude / Longitude with other data annotations as a Google Map
- Electrophoresis Gel Analysis
- Interpolation Calculator
- Photo Distance Object
- Auto resizing text annotations

Significant Defects & Issues Resolved

On Windows systems, Logger *Pro* files on read-only media (CDs etc) would fail to open due to a low level Microsoft defect. This has been corrected and Logger *Pro* now opens these files successfully.

Data synchronization with multiple sources is much improved and more consistent.

Logger Pro 3.4.2

International

Logger Pro 3.4.2 provides an upgrade path for 3.2.1 international users. International users should examine release note entries below for versions since 3.2.1 to gain an understanding of what has changed from 3.2.1 to 3.4.2.

Silent Installation

Logger Pro 3.4.2 provides and allows a silent install option for system maintainers of a large collection of machines for English only. For more information on the silent install option, please contact Vernier Software & Technology.

Issues Resolved

Over 80 issues were addressed by 3.4.2. The following list reflects some of the more notable issues resolved.

- Curve fits no longer fail when columns have single Greek character short column names.
- The correct curve fit now loads on opening all files saved with curve fits in versions before 3.4. This only affected saved curve fits after the polynomial fit in the curve fit menu.
- Removed a program failure that would occur if the user deleted the base column of a graph and stored the latest run.
- Changed the slope calculation to correctly account for blank cells.
- Text import now works correctly for utf8 text encoding (special characters).
- LabPro OS update no longer causes an intermittent failure.
- Curve fits with single-character x-axis labels now correctly use that label.
- Fixed program failure on keyboard entry when automatic linear fit helper object is the only object selected.
- Graph Options dialog now respects manual scale settings.
- Better LabPro driver install on Windows 2000.
- Fixed a LabPro power saving feature defect that prevented the unit from powering down in all situations.
- Fixed a Mac OS Logger Pro crash when a default folder name was left empty in the start-up file.
- Fixed an issue with FFT where removal of one input column accidentally removed the source.
- Fixed text edit issues in text boxes where some characters were turned into Greek characters.
- Fixed text select and replace cursor placement issues that made certain text manipulations difficult.
- Data Delete Column should now ONLY show columns that can be deleted.
- Corrected a situation that froze Logger Pro 3.4.1 using an analog sensor and a photogate sensor.
- Opening Digital Meter Options dialog no longer intermittently shows “invalid parameter” message.
- Fixed intermittent and rare file corruption problems.
- Fixed TI Graph Link silver cable problem when upgrading from Logger Pro 3.1.
- Non-numeric entries in numeric columns are no longer lost.

Logger Pro 3.4.1

Serial LabPro Communications

Fixed issue with slow (<0.2 Hz) data collection when using serial connection to LabPro.

Logger *Pro* 3.4

Supported Operating Systems

Support added for Mac OS X 10.4.1 and newer.

Signed Driver

LabPro USB driver for Windows XP is now signed by Microsoft. A signed driver can be installed without displaying a hardware wizard or requiring an administrator login. Initial Logger *Pro* installation still requires administrative rights, but connecting a LabPro to a new USB port will no longer require such rights.

Support for New Devices

- Support added for Go! Motion.
- Support added for TI-84 and TI-89 Titanium USB direct cables.
- Support for new sensors: Blood Pressure, Charge, Spirometer, and Hand Dynamometer.

User Parameters

User parameters allow control of values used in calculations, triggering, and analog output functions.

Control objects for user parameters have been added to the Insert menu. Parameter controls are page objects, much like a digital meter, that allow adjustment of parameter values using mouse or cursor keys.

Curve Fits and Modeling

Function Plotting. Modeling no longer requires the presence of data to draw a function on a graph, so a function may be plotted alone.

Modeling is improved. Manual curve fits may be modified after the fit is drawn on the main graph.

Curve fits now optionally stored as a calculated column of user parameters. Curve fits as parameters allows the fit to be manually adjusted after the fit.

Video Features

- Video capture from many USB camera devices, with automatic insertion of movies in Logger *Pro* pages, has been added to the Insert menu.
- When an experiment file containing inserted videos or pictures is created or saved, the referenced video or picture file is copied to the same file location.

- Video capture done with data collection will be inserted and already synchronized, ready for replay.
- Video analysis of a still photo.
- Rotated coordinate systems for video analysis allows one-dimensional analysis of motion that is neither vertical nor horizontal.
- Display of current frame time and frame number in videos.

Graph Features

Double Y graphs allow plotting of values with different magnitudes on a single graph, such as pH and temperature vs. time.

Adaptive point protectors for clearer printing of graphs with multiple traces. Point protectors are placed sparsely on runs as markers.

Better use of color for distinguishing runs. Optionally *Logger Pro* will use a single color in each run.

Time of day and Date in Time column and graphs has been added.

Optionally, horizontal axis labels may be rotated for graphs

There is a new option to make all rows major tics on graphs. For example, a data table listing twelve month names can be set to display all months on the horizontal axis.

Additional control of number displays added in data tables and graph labels. User can force scientific notation.

Styled Text

You can customize the style of the text (including adding Greek characters) as it appears in text boxes, annotations, and graph labels.

Data Collection

Automatic saving of backup data files in case of power loss or crash enabled.

Space bar toggles data collection on and off if the selected page object does not have space as a valid character. For example, if a graph is selected, pressing the space bar will start and stop collection, just like the Collect button. However, if a text object is selected pressing space will enter a space character. The space bar is a convenient way to start data collection.

New unit-changing paradigm allows changing units after data collection. Units are no longer tied to a particular calibration.

New Calculated Column Functions

- Event counting with photogates
- Modulo

- Boolean functions
- Collapse indirect for better display of photogate data
- Blood pressure calculations

Improvements

It is now possible to collect data more slowly than one point every 4.4 hours when LabPro is connected to the computer. Remote data collection continues to have the 4.4-hour upper limit to the time between points.

Improved speed for all Mac OS versions, especially OS 9. In particular, opening files with many calculated columns, or storing runs with many calculated columns, is much faster.

Improved flexibility in displaying photogate data without blank cells. The data table can hide rows based on any single column.

Improved examine mode for video analysis. Video now snaps to frame when a marked point is selected.

Additionally:

- Support added for QuickTime 7.
- Clarified sensor calibration dialog box.
- Improved printing of fine graph lines by setting gray grid lines to black.
- Sounds in videos are now enabled on all platforms.
- Improved tolerance of videos with uneven time steps, making more videos usable for analysis.
- The aspect ratio of an inserted video is preserved at file open.

Logger Pro 3.3

Support for Multiple Devices

Logger Pro 3.3 supports LabPro, Go!Temp, Go! Link, and Ohaus Scout Pro balances. Up to four devices may be used simultaneously, with no more than two LabPros. On Windows, only one Ohaus scale can be connected by USB. Note: USB devices connect automatically if detected. Serial LabPro (and Ohaus on Windows) require manual connection using Connect Interface from the Experiment menu.

Icons for On-Line Devices

The second row in the toolbar now holds icons for each on-line device, followed by any live readouts. These icons are buttons to summon the sensor setup dialogs for each device. Offline devices do not appear in the toolbar, but are available in the Experiment menu.

Keyboard triggering

Keyboard triggering is a new option that allows you to start data collection with minimal delay by pressing the space-bar or return key on the keyboard. Enable keyboard triggering on the Data Collection dialog, Triggering tab.

When keyboard triggering is active, clicking the Collect button merely prepares the system for data collection; collection begins at the next point when the space-bar is pressed.

Automatic Detection of Sensor Configuration

Logger *Pro* now detects the connection and removal of auto-ID sensors whenever live readouts are on and data collection is not active. When a sensor is added, a column is added to the Data Table. When a sensor is removed, its column is removed if the column is empty. To disable this behavior, turn off live readouts in the Experiment menu.

Automatic Updating of Page Objects

Logger *Pro* now automatically adds graphs or connects columns to existing graphs when sensors are added. This behavior simplifies setup, but it can also disrupt existing page layouts. As a result there are two new options to control automatic updates of page objects. You can separately enable or disable update page objects on File → New (the default is on) and after File → Open (the default is off).

Single-Point Calibration

Standard calibrations require two external references, but some sensors can be calibrated by setting an offset only, such as correcting a barometer for altitude. Logger *Pro* now allows single-point calibrations for any sensor that can be calibrated using two points.

Identify Column Data Source

Now that Logger *Pro* collects data from multiple devices it can be confusing as to where a given column is getting its values. You can now determine the source by device and channel number (e.g., LabPro 2, CH 3) by inspecting the Column Options dialog for a given column. In addition, the interface dialogs (Experiments → Set Up Sensors → Show all Interfaces) show the column name where the data from each channel is stored, in each sensor's Sensor Info display.

Finer Control of Live Readouts

All live readouts can still be disabled through the experiment menu but individual columns can be controlled through the column options. This can be very helpful to reduce the number of live readouts displayed on the tool bar. You can now turn on live readouts for the Motion Detector, although then the detector will click at all times.

Export to GIS Format

Logger *Pro* can export data in a form recognized by GIS software such as ArcView by ESRI.

New Page Objects

There are three new page objects. The Thermometer resembles a glass thermometer and can be used to display sensor readings. The Analog Gauge is a dial readout for sensor readings. The Animated Display is used to represent the motion of objects, usually from a Motion Detector column.

New Global Preferences

New preferences have been added. Optionally you can turn on a list of recently used files to be displayed in the File menu. This feature lets you quickly return to recent work, but some users may want to turn off the feature so that students cannot easily find other student work. The default is off. You can turn on large toolbar buttons so that *Logger Pro* will look nearly like *Logger Lite*. The default is set for standard *Logger Pro* buttons.

Supports New Sensors

The new Thermocouple (version shipping after May 2004) is now supported

Sensor Confirmation Dialog

The Sensor Confirmation dialog opens when you open an experiment file and *Logger Pro* does not detect all the necessary sensors. In version 3.2.1 and earlier we assumed that some non-auto ID sensors were present. Now we prompt, unless the Sensor Confirmation flag is cleared in the setting for that experiment file.

The Sensor Confirmation Dialog will disappear when you connect the requested auto ID sensors.

You can also indicate where any non-auto ID sensors are connected and manually close the dialog.

Logger Pro File Compatibility

Logger Pro saves files in the .cmbl file format. *Logger Pro* opens .xmbl, .gmbl, and .cmbl files as saved by *Logger Pro* version 3.0, 3.1, 3.2 and 3.2.1, and *Logger Lite*. *Logger Pro* will also open files saved in any version of Graphical Analysis 3 (.ga3 extension).

Other Changes

- Grouped graph behavior links axes and doesn't group the objects themselves.
- Page objects now click to the front like windows in Finder and Explorer.
- Sensor polling for new auto ID sensors is now on whenever live readouts are on.
- Option to hide helper objects added.
- Predictions now take on the units of the graph.
- If no interface is connected when *Logger Pro* is launched, there is no longer a prompt to connect an interface, since *Logger Pro* is often used for manually entered data.

Experiment Files

- Experiment files have been updated.
- Experiment files have been added for Workshop Physics and CPU (Constructing Physics Understanding).
- Experiment files for Workshop Physics, CPU, Interactive Lecture Demonstrations, Real Time Physics, and Tools for Scientific Thinking have been moved to a folder titled *Additional Physics*.

Logger Pro 3.2.1

Supported Operating Systems

Mac OSX 10.3 is now supported.

Enhancements

- Several printing issues corrected.
- Fixed problem with serial connection to LabPro and Events with Entry data collection.
- Two-Photogate timing mode corrected.
- Remote data collection for 12000 points and four sensors now works correctly.

Experiment Files

Experiment files have been updated.

International

Logger Pro 3.2.1 is more compatible with international options.

Logger Pro File Compatibility

- Logger Pro files created in version 3.2 and 3.2.1 cannot be opened in older versions of Logger Pro.
- Logger Pro files created in version 3.0 and 3.1 may be opened in Logger Pro 3.2 and 3.2.1.

Logger Pro 3.2

Video Analysis

Quantitative position information can be obtained from inserted videos. User can set scale, define an origin, and track the positions of up to three objects. Video analysis data may be synchronized with sensor data.

Strip Charts

A new strip chart graph type has been added. If the width of a strip chart is less than the data collection time, then the strip chart will scroll right to left as data collection or replay proceeds.

Rate Functions

The rate functions used for heart rate and respiration rate calculations has been improved to use a moving window, allowing more frequent updates to the reported values.

TI Connect Compatibility

Logger *Pro* 3.2 is compatible with TI Connect versions 1.1, 1.2 and 1.3, and will use the installed USB Graph Link driver for calculator data import, if one is present.

Optional LabPro OS update

You have the option to update the LabPro Operating System (OS) from 6.26 to 6.27 to obtain improved Photogate and remote data collection behavior. If your LabPro has an OS older than 6.26, Logger *Pro* 3.2 will require an update to 6.27.

Other changes

- New and better sensor select dialog, with groupings of sensors.
- Handles 60 samples a minute or hour and other unusual values with better display of time values.
- Comma decimal separator now supported in addition to North American standard period decimal separator.
- Additional Tool Tips added.
- Extra smart sensors are ignored when opening a file.
- Stainless Steel Temperature probe calibration now supported using three reference temperatures.
- Removed extra Motion Detector calibrations because you are able to remove automatically generated calculated columns now and they will not come back.
- Data can now start before the movie and still be synced.
- Improved Motion Detector calibrating.
- Large number support increased to 10^{50} .

Logger *Pro* 3.1

A major development in Logger *Pro* 3 is the addition of all the features of the Graphical Analysis 3 Program. Below is an overview of these features.

New Feature: Import From TI Device

You can use *Logger Pro 3* to transfer the data from a Texas Instruments graphing calculator to your computer for analysis or printing. Supported calculators are the TI-73, TI-73 Explorer, TI-82, TI-83, TI-83 Plus, TI-83 Plus Silver Edition, TI-84 Plus, TI-84 Plus Silver Edition, TI-85, TI-86, TI-89, TI-92, TI-92 Plus, and Voyage 200. To do this, you need a TI-GRAPH LINK cable, which is part of the TI-GRAPH LINK package and is sold by Vernier Software & Technology and other Texas Instruments dealers. This cable connects the TI graphing calculator to the serial or USB port of your computer.

New Feature: Import From Data Pro

You can import data from a Palm hand-held running DataPro into *Logger Pro 3* for analysis or printing by selecting Import from...DataPro from the File menu. The Palm Hot-sync feature will also send data to *Logger Pro*. You will need the DataPro package, which is sold by Vernier Software & Technology.

New Feature: Movies

You can synchronize data collection and a movie of the experiment so that as you replay the movie, the data are replayed on the graph. Supported movie formats are Video clip (.avi) and all file formats supported by the QuickTime player. For a complete list of supported formats, visit: apple.com/quicktime/products/qt/specifications.html

New Feature: Draw Prediction

Added to *Logger Pro 3* is the ability to predict what might happen in a given scenario. In certain lab experiments, you might want students to predict what is going to happen before they collect data. When you draw a prediction in *Logger Pro 3* on the graphs, the prediction will be saved in the data table for comparisons.

New Feature: Control Devices

Using *Logger Pro 3*, you can power devices with an output signal from the function generator built into LabPro. You can also power and control devices using the Digital Control Unit.

New Feature: Pages

Logger Pro 3 now has pages. This allows you to have different layouts of data in a single file. One page might contain a graph and table with five different data collections. Then on different pages you might perform comparisons between different collections or just go into each collection in more detail. Each page can almost act as a different file while still sharing data with all the other pages.

New Feature: Insert Additional Objects

Added to *Logger Pro 3* is the ability to add different objects besides the normal table and graph. These include rectangles and ovals. You can also add pictures to your experiment. Visual aids are especially useful for annotation, documentation, or instruction for a particular experiment.

New Feature: Data Browser

The Data Browser is a container for all the data in a *Logger Pro* file. In contrast, a given data table will not necessarily show all data columns.

For basic use of *Logger Pro*, you will not need to use the Data Browser. By default, all new data columns are added to all data tables. However, for more complex *Logger Pro* sessions, you may find it useful to display only a subset of the possible columns in a particular data table, or you may find it useful to not display a data table at all. In these cases, you will use the Data Browser to manipulate data columns and data sets.

What Changed from *Logger Pro 2* to *Logger Pro 3*

Below are questions you may ask if you are already familiar with older versions of *Logger Pro*. Also provided is a list of improvements and new features of *Logger Pro 3*.

Interfaces

Logger Pro 3 supports the Vernier LabPro interface, the Vernier Go! Temp, the Vernier Go! Link, and the Ohaus Scout Pro balance. You can import data from text files, Excel files, different TI Graphing Calculators using a TI-GRAPH LINK cable, a Palm OS hand held using *Data Pro* and the LabPro. *Logger Pro 3* will not work with the ULI, Serial Box Interface, or the MPLI.

LabPro and Collect buttons

The “Collect” button will no longer disappear when an interface is not connected. Instead there is constant feedback in the upper left-hand corner of the tool bar that tells the user whether or not the LabPro is connected. The Collect button will be faded if LabPro is not connected or recognized.

Windows

Logger Pro 3 now uses objects instead of windows. Every table, graph, meter, floating box and other window previously available in *Logger Pro 2* are now separate objects in *Logger Pro 3*. Each object can still be sized and moved. Using the Page menu items, objects can be aligned, layered, grouped, un-grouped, and locked on a page. To insert most objects, just select the object you wish to add from the Insert menu. Floating boxes, now called helper

objects, can still be found in the Analyze menu. Objects also have options, which can be opened by either double clicking on an object or selecting the desired object and choosing it from the Options menu.

Data Runs

Runs from Logger *Pro 2* are now called Data Sets in Logger *Pro 3*.

Calculated columns

Formula Columns are now called Calculated Columns in Logger *Pro 3*. Other than the name change, calculated columns act the same way as Formula Columns in Logger *Pro 2*, with some additional features like error bar calculations.

Supported File Formats

Logger *Pro 3* can open files in Logger *Pro 3* (.cmbl or .xmbl), Logger Lite (.gmbl) or Graphical Analysis 3 (.ga3) formats.

The Import Text feature, available from the File menu, can import data saved with the Export Text menu item or files created and saved from several other formats (e.g. Excel). The file must be in the tab-delimited text format used by Export Text (a .txt, .TEXT, .dat or .scv extension). Data are imported only into the latest data columns. You can also use this feature to import data prepared or collected in another program into Logger Pro.

Sensor Setup

The Setup menu has been removed. The Sensors Properties dialog has changed to the Sensor dialog, which can be accessed through the Experiment menu by choosing Set Up Sensors → Show All Interfaces from the Experiment menu. Each interface (LabPro, Go!Temp, Go! Link, or Ohaus balance) will have its own Sensor dialog.

If you are using an auto-ID sensor (such as a Stainless Steel Temperature Sensor or a Motion Detector), Logger *Pro* will identify your sensor at launch and load a default data collection mode. If your sensor does not auto-ID, choose an experiment file from the appropriate sensor folder or open the Sensor Setup dialog to manually configure Logger *Pro* for you sensor.

Mode and Sampling

In Logger *Pro 3* these two separate pages have been combined into the Data Collection dialog from the Experiment menu. As you change the mode for the collection, you will notice the options available will also change according to your choice. These options are similar to those in Logger *Pro 2*.

Repeat mode

In *Logger Pro 2* there was a separate Repeat mode for doing a Real Time Collection with repeats. In *Logger Pro 3*, Repeat is a checkbox option on the Time Based collection mode page of the Data Collection dialog.

Photogate Timing modes

Photogate Timing modes are not available in the Data Collection menu. Click on the Photogate Channel Box in the Sensor Setup window to change or select timing modes.

Logger Pro 3 also contains improved, more flexible functions (available in the New Calculated Column dialog from the Data Menu) to create calculated columns based on Photogate data.

Radiation Monitor setup

Once the Radiation Monitor is set up in the Sensor Setup window, the sensor status will be shown. The reading will be in Counts. The status of the Radiation Monitor is also displayed under the Toolbar. Click on the Radiation Monitor Channel Box in the Sensor Setup window to change or select timing intervals.

FFT graph and table

To insert an FFT graph, just select it from the Insert menu under Additional Graphs. In *Logger Pro 3* there is no longer a separate FFT Table; instead FFT data are placed in the regular data table. The FFT Graph Options dialog is similar to the FFT Options dialog found in *Logger Pro 2* except the two pages, Graph Features and FFT, have been combined into one page for easier access. FFTs have also been updated to allow more than one column to be plotted on a single FFT graph.

Histogram graph and table

To insert a Histogram graph, select Histogram from the Insert menu under Additional Graphs. In *Logger Pro 3* there is no longer a separate Histogram Table; instead histogram data are placed in the regular data table. The Histogram Graph Options dialog is similar to the Histogram Options dialog found in *Logger Pro 2*, but the Graph Features page has been renamed to Graph Options and the Axis Options page is now the Bin & Count Options page. Histograms have also been updated to allow more than one column to be plotted on a single Histogram graph.

Automatic/Manual Curve Fit

Both types of curve fits have been merged into a single dialog in *Logger Pro 3*. *Logger Pro 3* includes all *Logger Pro 2* curve fit options and many new features. For automatic curve fits, you are now able to define your own function for a curve fit. The preview window now acts like a normal graph with the ability to zoom and auto scale. To switch between a manual and automatic curve fit, just click the appropriate radio button located in the upper right of the dialog under Fit Type. After performing an automatic curve fit, you can switch to a manual curve fit to get an

even better fit if you desire. Sometimes it also helps to click the Try Fit button more than once to improve an automatic curve fit.

Auto scale axes

The functionality is now in the Axes Options tab of the Graph Options dialog, which can be called up by selecting Graph Options from the Options menu or by double clicking on the graph. Clicking on the axis labels on the graph may also change the axes. The scaling options for each axis are now independent of each other.

Printing

There is no longer an option to Print Screen. Instead, you can select Print from the File menu.

Remote data retrieval

Remote data retrieval from the LabPro is now in the Experiment menu. The functionality remains the same except that when you setup the LabPro for remote data, you can now print out your collection settings. Retrieval of data from the LabPro remains the same.

Calibration Folder

The calibration for sensors has been updated and is now stored in a database, so it's no longer necessary to supply a calibration folder.

Experiment Folder

The experiment folder is now located within the default folder, which can be defined in the Preferences dialog.

Over Range Auto-scale option

The Over Range Auto-scale option in the Graph behavior preferences page in the Preferences dialog has been removed. In Logger *Pro* 3 you can get the same effect by going to the Axes Options, which can be accessed either through the Graph Options (from the Options menu or double clicking on the graph). There are scaling options for each axis, so the auto-scaling can be independent of each other.

Default Interface

Logger *Pro* 3 allows the creation of offline interfaces, replacing the single default interface of Logger *Pro* 2.

Options menu

The Options item found in the Experiment menu in *Logger Pro 2* has moved to the File menu and renamed *Settings for "this file."* You will find all of the experiment options in this dialog except the Enable automatic Curve Fit check-box, which has been moved to the Preferences dialog. This new Settings dialog will save your options for the entire file, even after closing the file.

Store Latest Run

To store a set of data, choose Store Latest from the Experiment Menu.

Information about a file

The *About "this file"* item in the Help menu has been removed. In *Logger Pro 3* you can place information about a page in the Page Options item in the Page menu. You can set the Page Information text in the Page Options dialog for page 1. You can also choose for it to be displayed when the page is first viewed, or when the file is opened.

Undo function

In *Logger Pro 2* there was an undo command in the Edit menu. *Logger Pro 3* expands on this command and allows you to undo almost any action that you perform. Collecting data, deleting or inserting an object, or changing the scale of the graph, for example can all be undone. If you perform too many undo operations, you now have the option to perform a redo that essentially undoes an undo.

Text Annotation

Text Annotation in *Logger Pro 3* has been updated with many new features. For one, each text annotation can be connected to a point, thus describing information about what is happening to the graph at a particular point in time. As the graph is moved, the text annotation will stay with its point. You can also have more than one text annotation for a single point.

Stretchy Axes

Just outside the X or Y-axis line, you can use the cursor to manually rescale or stretch the axis.

LabPro Operating System Changes 6.26 (required) and 6.27 (optional)

Logger Pro 3 requires at least LabPro operating system 6.26. If a LabPro is connected with an earlier version, a dialog box will open offering to update the LabPro. You must update the LabPro to use it with *Logger Pro 3*. Beginning in version 3.2.1, *Logger Pro* will update the LabPro OS to version 6.27. This update is optional, and be applied by connecting a LabPro, clicking the LabPro icon in the toolbar, and clicking the Update LabPro OS button in the Sensors dialog for that LabPro.

In version 6.26 sensors were allowed to warm up after pressing the START/STOP button. In 6.27 the sensors are warmed up at setup time. After any warmup period has elapsed, data collection begins immediately on pressing START/STOP.

Version 6.27 updates and fixes:

- The build has DataMate 6.15 that includes the TI-73 Application.
- Fixed problem where communications were ignored during sensor warm up period.
- Fixed problem of photogate timebase and analog timebase being out of sync.
- Fixed problem of a first point of monitored (not collected) data being max value when using high speed data collection (DataMate),
- Fixed problem of corrupted data if photogate state transitions during start of pendulum experiment
- Fixed problem of using two Radiation monitors at the same time (*Logger Pro*)

LabQuest Operating System Changes 1.0 (required) and 1.1.1 and newer (recommended)

Logger Pro 3 requires at least LabQuest operating system 1.0 however it is highly recommended that you upgrade to 1.1 or newer. If a LabQuest with 1.0 software is detected, a warning will be displayed advising that the LabQuest unit be upgraded. This message will be displayed every time the LabQuest is detected, but the LabQuest can still be used for most applications. If a Logger Pro feature requiring newer LabQuest software is used, unexpected results such as freezing, crashing, or data loss may occur. We strongly recommend updating LabQuest to at least 1.1.1 (for 3.6 international users) or 1.3 for all others or newer. See www.vernier.com/labquest/updates/ for information on updates.

Version 1.1 Notes:

- Rotary Motion sensor and Radiation Monitoring are supported.
- Internal Microphone is supported.

Version 1.1.1 Notes:

- This update corrected an issue that could lead to the LabQuest “freezing” for some hardware units.

Version 1.2 Notes:

- GPS support.
- 100khz data collection
- File management and interchange improvements
- This LabQuest may have difficulty opening QMBL files created by *Logger Pro* 3.6i (international version). 3.6i users may wish to work around this issue by changing the version tag information in the QMBL from “3.6i” to simply “3.6” or “3.6 i” (with a space). Alternatively, 3.6i users can upgrade to 3.7i (when available) where this issue is fixed. Only users attempting to use QMBLs **created by *Logger Pro* 3.6i** have this issue.

Version 1.3 Notes:

- Supports the GC-Mini
- Supports Watts Up?
- Supports Power Amp and Function Generation
- This LabQuest may have difficulty opening QMBL files created by Logger *Pro* 3.6i (international version). 3.6i users may wish to work around this issue by changing the version tag information in the QMBL from “3.6i” to simply “3.6” or “3.6 i” (with a space). Alternatively, 3.6i users can upgrade to 3.7i (when available) where this issue is fixed. Only users attempting to use QMBLs **created by Logger Pro 3.6i** have this issue. Version 1.4 Notes:
-

Version 1.5 Notes:

- Supports data marking
- Supports the anemometer, Melt Station, and 30V sensors
- Enhanced battery power management.
- Enhanced networking.
- Data Matrix data collection mode
- WDSS support

Vernier Software & Technology
info@vernier.com
www.vernier.com